

Prepable Titanium Abutments

Prepable Titanium Abutments are available in various angulations, heights and range of anatomical shapes for all positions and match perfectly the Tissue Formers and the natural profile developed during healing. Esthetic abutment designs minimise the need for further adjustments and allow an easy, fast and esthetic restoration.



- *Intended use* – Cement-retained single and multi-unit restorations.
- *Material* – Commercially Pure Titanium Grade 4.
- *Packaging* – Non-sterile; includes abutment, laboratory screw and abutment screw.
- *Esthetiline Solution* – Prepable Titanium Abutments are part of the Esthetiline solution. The shapes match the profile of the Tissue Formers and an optimal result is achieved by choosing the Prepable Abutment that matches the Tissue Former placed at surgery.
- *Abutment modification* – the abutments may be modified chairside, intra- or extra- orally, or by the technician on a laboratory model using Titanium Abutment Preparation Kit, tungsten carbide or diamond burs.
- *Indexing* – It is recommended that the implant is indexed prior to placing the non-rotational abutment to ensure optimal placement and minimal preparation.
- *Procedure* – The abutment is screwed into the implant and the restoration is placed over the prepared/modified abutment. For more information please refer to the Neoss Implant System Guidelines.

Assortment

Item no.	Description
31280	Prepable Ti Abutment 1.5mm 0° - Wide Incisor
31281	Prepable Ti Abutment 1.5mm 15° - Wide Incisor
31282	Prepable Ti Abutment 1mm 0° - Narrow Incisor
31283	Prepable Ti Abutment 1mm 15° - Narrow Incisor
31284	Prepable Ti Abutment 3mm 0° - Narrow Incisor
31285	Prepable Ti Abutment 3mm 15° - Narrow Incisor
31286	Prepable Ti Abutment 3mm 0° - Canine
31287	Prepable Ti Abutment 3mm 0° - Pre-molar
31288	Prepable Ti Abutment 3mm 0° - Molar
90136	Titanium Abutment Preparation Kit

Note: 31282 - 31285 correspond to previous articles 31208 - 31211

